STERILCONTAINER™ System

US

INTRODUCTION:

The Aesculap STERILCONTAINER System is a reusable rigid container system used for the packaging, transportation, and storage of instruments prior to, during, and after sterilization. It consists of various sizes of lids and bottoms, with assorted accessories such as baskets, filters, indicator cards and tamper proof locks.

The Aesculap STERILCONTAINER System is an alternative to the traditional woven or nonwoven wrapping materials used to package surgical instruments and other supplies for sterilization.

CARE & HANDLING:

The Aesculap STERILCONTAINER System is manufactured of specially formulated aluminum alloy. The surface is protected by a layer of anodized oxide to prevent corrosion. If properly cared for using the instructions in this booklet, your Aesculap STERILCONTAINER can provide many years of instrument protection and sterility assurance.

DO NOT USE ABRASIVE CLEANERS, METAL BRUSHES OR ABRASIVE CLEANING PADS. USE OF ABRASIVE PRODUCTS CAN CAUSE PERMANENT DAMAGE TO STERILCONTAINER SURFACES. USE OF ABRASIVE CLEANERS OR PADS WILL RESULT IN WARRANTY EXCLUSION.

CLEANING:

Recommended Cleaners:

Use only mild alkaline, sodium carbonate-free, neutral pH (7) detergents to clean effectively without causing damage to the anodized layer of aluminum. If in doubt, contact the detergent manufacturer or supplier to determine suitability for use in cleaning the Aesculap STERILCONTAINER.

Pre-Cleaning Preparation:

1. Remove the lid from the container bottom.
2. Remove the basket and any instruments from the container.
3. Remove the lid retention plate(s).
4. Remove the container bottom retention plate(s) if using perforated container bottoms.
5. Discard the disposable filter(s).
6. Remove all processing indicators and disposable locks.

NOTE: Aesculap disposable non-woven filters are designed for one-time use only and should be discarded after each processing cycle.

Cleaning:

Manual Cleaning:

1. Use a soft sponge and a mild detergent and clean the STERILCONTAINER and all the components under water.
2. Rinse thoroughly under running water to remove all detergent residue, as residual detergents can affect the anodized oxide layer of aluminum.
3. To remove sterilization adhesive tape remnant of surface abrasions, we recommend the use of Aesculap-Eloxal Cleaner (Catalog number JG601). This is a non-abrasive cleaner. Apply the cream with a soft dry cloth and rub to polish the surface. Thoroughly rinse the STERILCONTAINER under running water to remove all residual cleaning cream.
4. Thoroughly dry the STERILCONTAINER and all components with a soft dry cloth.

Note: *Wear proper protective personal attire when cleaning the STERILCONTAINER.*

Mechanical Cleaning:

1. Place the STERILCONTAINER bottom in the washer with the inside surface facing down to avoid water collection.
2. Fold the handles towards the inside of the lid. Place the lid with the inside surface facing down to avoid water collection.
3. Retention plates should be placed away from the direct force of pressurized washer jets to avoid damage during the washing cycle. The wire baskets provided by the mechanical washer manufacturer should be used to hold STERILCONTAINER components during the automated cleaning cycles.
4. Thoroughly dry (either with a soft, dry cloth or air dry) the STERILCONTAINER and components before final assembly.

If you have further questions regarding cleaning practices for the STERILCONTAINER System, contact your local Aesculap Sales Representative or the Customer Service department at 1-800-282-9000.

ASSEMBLY FOR USE

Pre-Assembly Inspection and Preparation:

1. Identify surgical instrument set for sterilization.
2. Identify appropriate size basket and STERILCONTAINER, as well as STERILCONTAINER bottom type (perforated or solid) for the method of sterilization chosen. **(It is not recommended to use solid bottom containers in gravity steam sterilizers.)**
3. Assure STERILCONTAINER component pieces are completely dry.
4. Inspect the rim of the lid to ensure the gasket is in good condition and free from cracks. A cracked gasket indicates age and/or deterioration and should not be used. Remove the lid from service and return for repair.

NOTE: An optional inner drape may be used to wrap the basket and basket contents to facilitate the formation of a sterile field when opened.

Filter Assembly:

1. Place one sheet of the appropriate Aesculap filter paper over each perforated section on the inside of the STERILCONTAINER lid. When using Aesculap filter paper with indicator ink (US751), be sure indicator dot faces toward the inside of the container when the lid is in place and the indicator dot is visible through the retention plate.
2. Secure each filter with the retention plate designed for use with the STERILCONTAINER lid.
3. Perforated Bottoms: Place one sheet of the appropriate Aesculap filter paper over the perforated section in the STERILCONTAINER bottom. Secure the filter with the retention plate designed for the STERILCONTAINER bottom.
Solid Bottoms: This model has no perforated section and does not require a filter or retention plate.

Instrument and STERILCONTAINER Assembly:

1. Sort and assemble thoroughly cleaned and dried instruments into the instrument basket(s), according to established hospital procedures.
2. Place assembled instrument basket(s) into the prepared STERILCONTAINER bottom.
3. Place assembled lid onto the STERILCONTAINER bottom, aligning handles on bottom with latches on lid.
4. Simultaneously close both locking latches on the STERILCONTAINER lid.

NOTE: Aesculap has validated a weight of 16-20 pounds for baskets and basket contents. Hospitals should refer to AAMI/AORN guidelines on weights and weight limit. All instruments should be assembled to allow for uniform exposure to sterilization agents.

CAUTION:

Leave one inch of free space between the instruments and the rim of the container bottom for effective processing. Basket handles may encroach into this clearance space as long as they do not interfere with the lid's filter retention plate or lid closure.

Processing Assembly:

1. Select the appropriate Aesculap Indicator Card and insert into the holding bracket on the outside of the STERILCONTAINER. A tab at one end of the indicator card will facilitate insertion and removal.
2. Insert the appropriate Tamper Proof Seal into the locking channel on each end.
3. Secure and lock the seals.

NOTE: Use of internal and external indicators should be in accordance with in-house protocol, determined by the user.

Loading the Sterilizer:

1. **In all methods of sterilization**, the STERILCONTAINER should be placed **flat** for effective sterilization and drying.
2. **In all methods of sterilization**, the STERILCONTAINER should be positioned on the autoclave cart below wrapped items for optimum sterilization and drying conditions.
3. **Stacking:** Recommended for **high-vacuum cycles only**. Stacking should not exceed 16-18" in height for effective air removal and adequate steam penetration. Both solid and perforated bottom containers can be stacked.

PROCESSING:

1. Run loaded sterilizer according to time and temperature recommended by the sterilizer manufacturer for chosen cycle.
2. At completion of the cycle, the sterilizer door should be opened approximately 6" and left open for a period of 15 minutes to aid in drying, as recommended by the sterilizer manufacturer.
3. The STERILCONTAINER needs to remain on the container cart, in a draft free area, until cool enough to handle. Cooling in a draft-free area minimizes the formation of condensation on the inside of containers. If condensation on the inside of containers comes into contact with outside contaminants, sterility may be compromised. Furthermore, condensation inside containers could cause rust to form on some instruments.

NOTE: A STERILCONTAINER with a solid bottom may require additional cooling time; the additional required cooling time to be determined by the user.

SUGGESTED STERILIZER CYCLE PARAMETERS:

The following parameters are based on the generally accepted recommendations of most sterilizer manufacturers. Each facility may need to run internal testing to determine if adjustments are necessary for their facility.

STEAM CYCLES:

Gravity:

CAUTION: DO NOT USE SOLID BOTTOM CONTAINERS IN GRAVITY STEAM CYCLES.

Temp:	250° F
Exposure Time:	30-60 Minutes (depending on load contents)
Cycle Dry Time:	15 Minutes (minimum)
Cool Time:	Varies according to load contents
CAUTION:	Cool drafts from air ducts or other air currents should be avoided during the cooling phase to avoid post-sterilization moisture caused by Rapid Cooling Syndrome.

High Vacuum (pre-vacuum, three pulse, standard):

Temp:	270° F
Exposure Time:	4 Minutes (minimum)
Cycle Dry Time:	15 Minutes (minimum)
Cool Time:	Varies according to load contents
CAUTION:	Cool drafts from air ducts or other air currents should be avoided during the cooling phase to avoid post-sterilization moisture caused by rapid cooling syndrome.

ETHYLENE OXIDE CYCLE:

Temp:	125° F - 130° F
Pre-Vacuum:	25" Hg (minimum)
Humidity:	40-60% RH
Gas Pressure:	600 mg/L (minimum) EO Gas Mixtures may vary ie. 10/90%by weight or 100%
Exposure Time:	105 Minutes (minimum)
Post Vacuum:	20" Hg (minimum)
Aeration:	8 Hours (minimum)
CAUTION:	Articles contained within the STERILCONTAINER may require longer aeration periods; adequate aeration times to be determined by the user.

NOTE: Consult with the manufacturer of the sterilizer for specific recommendations.

THE STERILCONTAINERTM SYSTEM WARRANTY:

The Warranty covers the following conditions:

The Aesculap STERILCONTAINERTM System is guaranteed to be free of functional defects in workmanship and materials when used appropriately for its intended purpose. Any STERILCONTAINER delivered from Aesculap, Inc. proving to be defective will be repaired or replaced, at Aesculap's discretion, at no charge to the customer.

Periodic Maintenance:

Periodic Maintenance is suggested to ensure all parts are in good working order. Proper care will enhance the lifetime of your STERILCONTAINER System.

Repair or Replacement:

Under this Warranty, and at its sole discretion, Aesculap will repair or replace any parts of the product found to be defective in materials or workmanship.

Aesculap will determine, at its sole discretion, the cause and nature of any defect, the necessity and manner of repair or replacement, and all other matters pertaining to the conditions of the STERILCONTAINER System.

Aesculap, Inc.'s Returned Goods Policy applies to all returned items (refer to details in current price list).

Exclusions:

This Warranty shall not apply to any conditions(s) or damage resulting from negligence, misuse, improper handling, improper cleaning, improper opening techniques, or unauthorized repair work including but not limited to:

- Incidents of abuse such as denting of the STERILCONTAINER due to dropping or other instances of mechanically applied pressure.
- Damage determined to be related to caustic or abrasive cleaning agents.
- Color-fading due to harsh detergents or abrasive cleaning agents or materials.
- Items modified by the customer.
- Items modified or customized by Aesculap at the request of the customer.
- Damage from fire, flood and other occurrences not under the control of Aesculap, Inc..

The Warranty period for the STERILCONTAINER gasket is three full years from the date of purchase.

Questions:

Please direct your questions concerning this device to:

Aesculap, Inc.

3773 Corporate Parkway
Center Valley, PA 18034

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- Daños derivados de incendios, inundaciones y otros acontecimientos que escapan al control de Aesculap, Inc..